

## 1200 Stretcher Lift

---



The 1200 Stretcher Lift indexes seamlessly with the 6302G series of supine tubs for secure comfort.



An optional digital scale (Model IH1100) is accurate to nearest pound or kilogram.

### **Full-body support for a fully relaxing bathing experience.**

The 1200 Stretcher Lift is designed for use especially with the 6302G series of supine tubs. For residents requiring assistance, this lift features stretcher wings that can be raised or lowered from either end for maximum flexibility in transfer and bathing of residents.

The mesh stretcher cover, fitted to a stainless steel frame, provides comfortable and contoured full-body support. As a result, residents can relax and enjoy the bath while caregivers can rest assured that residents are secure and safe during the bath.

Constructed of durable heavy-gauge steel and featuring hospital-grade casters, the 1200 Stretcher Lift features a vertical range of 30 inches and a weight capacity of 350 pounds to accommodate a wide range of residents. A precise digital scale is available as an optional accessory.

