



2600KL

Wheelchair ramp scale



2610KL

Wheelchair scale with dual ramps



Advantages

- ✓ **EMR Connectivity**
Scale interfaces with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's Body Mass Index (BMI)
- ✓ **Versatile Platform**
Wide and low profile platform accommodates wheelchairs and also provides a comfortable weighing experience for bariatric patients
- ✓ **Interactive Display**
High-contrast color screen provides a user-friendly interface
- ✓ **Casters and Folding Ramps for Easy Portability and Storage**
4 wheels including 2 casters afford ease of movement and the folding ramps make it simple to store
- ✓ **On-Screen Help Menu**
Guides users through operating functions and provides easy access to multiple settings
- ✓ **Everlock® (patent pending)**
Option that allows users to permanently lock the scale in either LB or KG to improve safety
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost



Specifications

- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x h x d): 36" x 32 1/4" x 2 1/2" / 914 mm x 819 mm x 64 mm
- Ramp Size (w x d): 32 1/4" x 9 7/8" / 819 mm x 248 mm
- Product Footprint 2600KL/2600KG & 2600KL-BT/2600KG-BT (w x h x d): 48 1/8" x 43 7/8" x 49 3/8" / 1222 mm x 1114 mm x 1254 mm
- Product Footprint 2610KL/2610KG & 2610KL-BT/2610KG-BT (w x h x d): 56 3/4" x 43 7/8" x 49 3/8" / 1441 mm x 1114 mm x 1254 mm
- Product Weight 2600KL/2600KG & 2600KL-BT/2600KG-BT: 93 lb / 42 kg
- Product Weight 2610KL/2610KG & 2610KL-BT/2610KG-BT: 98 lb / 44 kg
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Wheels: 4 (2 casters)
- Power Source: 100-240V adapter included, 6 D batteries (not included)
- Connectivity: USB, optional Pelstar wireless technology
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Everlock®
- 2 year limited warranty

Model	Description
2600KL	Digital wheelchair ramp scale, pounds & kilograms
2600KL-BT	2600KL with Pelstar wireless technology
2600KG	Digital wheelchair ramp scale, kilograms only
2600KG-BT	2600KG with Pelstar wireless technology
2610KL	Digital wheelchair with dual ramps, pounds & kilograms
2610KL-BT	2610KL with Pelstar wireless technology
2610KG	Digital wheelchair with dual ramps, kilograms only
2610KG-BT	2610KG with Pelstar wireless technology
SS-2600KL	ScaleSurance 2 year extended warranty for 2600KL, 2600KG, 2600KL-BT & 2600KG-BT*
SS-2610KL	ScaleSurance 2 year extended warranty for 2610KL, 2610KG, 2610KL-BT, & 2610KG-BT*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
B2600RAMP	Optional secondary ramp for 2600KL, 2600KG, 2600KL-BT, & 2600KG-BT
PELSTARONE-P	Wireless thermal printer