

## TheraPure™ 6302G & 6302GADH Supine Tubs

---



An integrated console provides the operator with easy access to the auto-fill/auto shut-off and anti-scald mixing valve, thermometer and whirlpool jet controls.



Adjustable height model (ADH) reduces strain on caregiver by raising or lowering tub at the touch of a button.

### **Pipeless™ technology meets flawless style.**

The TheraPure 6302G and 6302GADH supine tubs combines the purity and benefits of pipeless technology with updated contemporary styling, including a smooth fiberglass underside.

Large enough to allow full immersion for resident comfort, these tubs are designed to maximize caregiver ergonomics for ease of use and safety.

Additional space in the shoulder and hip areas make resident turning easier, while the narrower side walls allow the caregiver to be closer to the bather. The tubs include an auto fill feature and adjustable temperature settings. The adjustable height model 6302GADH can be raised to a comfortable working height for the caregiver, reducing back strain.

Four whirlpool jets provide relaxing hydromassage for semi- and non-ambulatory residents. The Invacare 1200 Stretcher Lift is designed to work in harmony with the 6302G series tubs.



