

EMR



With the highest capacity of Health o meter® Professional's remote display scales, the 752KL is a great supplement to accommodate ambulatory bariatric patients. The remote display and portable base make it convenient to move the scale anywhere it is needed.

Advantages

- ✓ **EMR Capable**
 Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions
- ✓ **Remote Display**
 The remote display provides versatility with the ability to be mounted on a wall or placed on a table
- ✓ **600 lb Capacity Floor Scale**
 The 752KL is the market's highest capacity floor scale, providing more versatile use
- ✓ **Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **Calculates Body Mass Index (BMI)**
 Scale functions include the ability to calculate the patient's BMI



Optional Carrying Case
(order # 64771)

Remote Display can be
mounted on a wall or
placed on a table



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
14 1/4" x 14 1/4" x 2 5/8"
362 mm x 362 mm x 67 mm
- Display: 1" / 25 mm LCD
- Connectivity: USB, optional Bluetooth[®]
- Power Source: 100-240V adapter included,
6 C-cell batteries (not included)
- Optional Carrying Case: 64771
- Functions: LB / KG Conversion, LB / KG Lock Out,
BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- 2 year limited warranty

Model	Description
SS-752KL	ScaleSurance 2 year extended warranty
ADPT50	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
64771	Optional carrying case for 752KL
BLUETOOTHONE-P	Wireless thermal printer