



The Astor-Bannerman Nivano range of wall mounted height adjustable changing tables have been a popular choice among families and Occupational Therapists for many years. This is because they are easy to use, comfortable, reliable and look great in any room without having a clinical feel.

However, the new improved Nivano 2 adjustable shower/ changing bench offers even more benefits. It is simple to operate from a handset control and comes in a range of sizes from 1000mm right up to 1900mm to suit both paediatric and adult applications. This gives the widest choice of stretcher lengths on the market today - at no extra cost.



Shown with optional guard

With a weight capacity of 150kg (220kg model available), the Nivano 2 has a comprehensive range of options so it can be tailored to suit your needs.

The durable stretcher material is comfortable as well as quick and easy to clean, saving precious time. This, along with part stainless steel construction means the Nivano 2 will give many years of worry free use, even in wet and humid or multi-user environments, especially when combined with ongoing maintenance for extra peace of mind.

The 600mm travel range enables the user to be gently lifted to a safe and comfortable working height for washing, showering, changing and dressing. The platform stows neatly when not in use, held in place with twin gas struts, and features battery back-up for safe lowering in the event of a power failure.

The Nivano 2 can be tailored for you with the addition of optional extras adding comfort, convenience and unique customisation.



With optional tilting head/backrest



Optional guard with padded bumper



Stows neatly against the unit when not in use



With optional removable mattress



STANDARD FEATURES

- SWL 150kg/24st
- Height range 300mm - 900mm (600mm)
- Height adjustable (mains electric)
- Handset operation
- Battery backup
- Floor and ceiling hoist access
- Part stainless steel construction
- Durable non-rip high tension vinyl stretcher
- Suitable for use in shower areas and wet rooms
- Protrudes just 242mm from wall when stowed
- 3 Year warranty*



OPTIONAL EXTRAS

- Stretcher sizes from 1000-1900 x 700mm
- SWL 220kg/35st available (1700-1900 models)
- Battery operated model (CTB2)
- Fold away safety guard
- Protective padded guard bumper
- Anti-splash skirt
- Tilting head/backrest (1700-1900 models)
- Protective foam wedge support
- Removable mattress
- Special fully coloured models available
- Range of service & maintenance packages - *Contact for more information*

*Terms and conditions apply






Model: Nivano 2 (CTE2) Electric Height Adjustable Changing Table

Product Code/s: 00.038.11.1000 - 00.038.11.1900 (150kg)
00.038.12.1700 - 00.038.12.1900 (220kg)

Stretcher Size: 1000mm - 1900mm x 700mm (150kg)
1700mm - 1900mm x 700mm (220kg)

Height Range: 300mm - 900mm x 700mm (600mm)

Features:

				
Height Adjustable	Mains Electric	Battery Backup	150kg/24st SWL	3 Year Warranty

Options:

					
Battery Model	1000 - 1900mm	220kg/35st SWL Model	Tilting Head/Backrest	Service Packages	Many More Options

Dependency Level:

High Dependency

Application:

Private Home	School	Leisure Centre	Care Home	Nursing Home	Respite Home	Hospital	Changing Places 
--------------	--------	----------------	-----------	--------------	--------------	----------	---

Technical Drawing:

