



Yes, you can.®

Invacare® Carroll CS Series™ CS7™ Bed

Bed featured with ThinkSoft® Positioning Device, Drake bed ends in Solar Oak and Invacare® Solace® mattress



The innovative CS7 Bed has become the industry standard for long-term care beds. With an Auto Contour feature that simultaneously raises both the head and knees, residents stay comfortably in place while being better positioned for reading, watching TV or talking to family members. The simultaneous movement also helps prevent painful skin shearing that might otherwise occur during repositioning.

Invacare® Carroll CS Series™ CS7™ Bed



Tool-less Assembly
Fast connect system allows for faster set-up and easier alterations.



Hand Control Pendant
A rugged, water-resistant pendant with fewer recessed areas is easy to clean and very intuitive for the resident to use.



Mobility On Demand
Floor lock system secures the bed to the floor at all height ranges but allows for mobility on demand.



Slat Deck
Exclusive arched slat deck design allows for improved infection control and reduces wear and tear on mattresses.

CS7 Bed Specifications

- 7" – 30" Height Range
- 500 lbs. Weight Capacity
- Auto Contour
- One-Touch Floor Lock
- I.V. Rod Holders – standard
- Meets UL60601-1
- Meets CAN/CSA 22.2 No. 601.1-M-90
- Meets CAN/CSA 22.2 No. 601.1-2-38:03
- 8 Button Hand Pendant
- Ergonomically Designed, Rugged, Water-Resistant Pendant
- Universal Deck Length (76" or 80")
- Mattress Keepers
- Optional Attendant Control Panel
- Optional Battery Backup
- Ratchet Heel Lift

CS7 Bed Accessories

- Trapeze
- Hi-Low Lockout
- CPR Release
- 4" Bed Extender
- Safety Package

Invacare Corporation

www.invacare-cc.com

info_icc@invacare.com

USA
One Invacare Way
Elyria, Ohio 44035
800-668-2337

Canada
994 Hargrieve Road
London, Ontario
N6E 1P5
800-668-2337

Trademarks are identified by the symbol™ and®.
All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

© 2011 Invacare Corporation. All rights reserved.
Form 08.30.9/11



Yes, you can.®