

Aluminium sling lift CARLO

The versatile patient lift for all areas of care ...

BEKA Hospitec develops and produces high quality and user friendly sling lifts, which offer your resident a maximum of safety and comfort and at the same time provide for optimal working conditions for your staff ...



Aluminium sling lift CARLO Comfort EP

The transfer lift offers the highest comfort for the work of the nursing staff. Due to the light aluminium construction, the lift is extra flexible and easy to move. The balance bow enables a safe and comfortable transfer of the resident. By means of the electrical positioning system the resident can be effortlessly moved into a sitting as well as into a lying position. No unsettling swinging movements or vibrations occur during lifting or transfer. This lift is designed for the application of 4-point-slings with safety clip system.

The CARLO Comfort EP is characterized especially by its optimal and safe weight balancing with a maximum of operational comfort.

Product advantages CARLO Aluminium sling lift

- 30% lighter construction due to 100% aluminium parts (carrier frame, lifting column, sling bow)
- Applicable for the shower
- Extra flexible and easy to move
- 4 double castor, 2 of them lockable
- Electrical height adjustment with lowering of the belt bow to the floor
- Electrically spreadable carrier frame
- Hand control with display of charging condition of the battery
- Locking device for fixation of balance-bow in moving direction
- Adapter for insertion of different bows
- Manual emergency lowering system
- Max. patient weight 185 kg
- Optional: calibrated patient scale

The optimal solution

for your nursing staff and your residents ...



The CARLO patient lifts are the ideal devices for lifting and transfer of persons, who need a high level of care, who can not support themselves or are altogether limited in their mobility ...



Aluminium sling lift CARLO Comfort GP

With the balance bow, the resident can be turned during transfer in an ergonomic and comfortable position, with the effect of an optimal pressure relief and weight balancing. The adjustment of the bow is facilitated by a gas spring and offers the possibility to stop in every desired position. This lift is designed for the application of 4-point-slings with safety clip system.

The 100% aluminium construction offers a simple and easy manoeuvrability.



Aluminium sling lift CARL Classic

The multi-purpose transfer lift for a safe and comfortable transfer of persons limited in their mobility. The adaption bow is designed for sling-belts. The basic frame, lifting arm and sling suspension are made of aluminium.

