

Lifting tub AVERO Classic

The classic, elegant nursing tub ...

The AVERO Classic lifting tub meets the highest demands on quality and reliability. The application of appropriate seat- and stretcher-lifts is possible without any problems.



Lifting tub AVERO Classic

The specific characteristic of the AVERO Classic tub is its timeless design, which is reflected in a pure, simple outside panelling and in a very beautiful interior tub space in body shape. An application of suitable seat- and stretcher-lifts is possible without problems, whereas special attention was paid to the fact that the resident has sufficient water volume for a relaxing bath.

The AVERO Classic provides unrestricted access for all hoists, therefore, even in small rooms and narrow spaces one can use bathing lifts without problems. The adapted back shape of the tub body with the distinctive arm rests and the premium, chromed brass grips provide a balanced tub interior shape, which allows the resident a relaxed stretch out in the tub.

Product advantages lifting tub AVERO Classic

- Easy operation by symbols
- Stepless height adjustment for optimal working height
- Pore-free surface for easy disinfection and optimal hygiene
- Unrestricted access for all hoists
- Ergonomical interior tub body
- Arm rests with grips for more comfort and safety
- Domestic appearance
- Scalding protection for hand shower and tub spout
- Optional: air-jet-whirl-system, disinfection system, automatic tub filling stop, colour-light-therapy-system, music-sound-system, tub shortener, aroma-therapy-system, colour at choice