

Invacare Roze Stand-Up Lift

Premium power stand-up lift for patients up to 450 lb.

Product ID: ROZE

Product Description

The Invacare Roze Stand-Up Lift is designed to promote independence for weight-bearing users up to 450 lb. as well as those needing rehabilitation support. Easy to use, flexible and versatile, this battery-powered lift leaves lifting to the lift for caregiver safety.

Features

- Height-adjustable knee pad for user comfort
- Large push handle and simple sling attachment simplify everyday transfers
- Removable footplate offers flexibility for users only needing help rising
- Small footprint combined with lightweight construction provides easy maneuverability
- Anti-entrapment motor safety device

Specifications

Height Adjustment:	Height At Sling Hook Up: 37"-63.5" H
Dimensions:	Base Dimensions: 4.8" H Clearance 44" W Open 26" W Closed
Caster Options:	4"
Product Weight:	117 lb.
Product Weight Capacity:	450 lb.
Battery:	Battery Life Indicator: Audio, Visual Battery: 24V Rechargeable 2w/Unit
Safety System:	Emergency Controls: Stop, Up and Down
Charger:	Charge Time: 6 hrs. Charger Input: 100-240V AC Charger Output: 29.5 DC Max
Warranty:	Limited Warranty - 5 Year - Lift; 2 Year - Electronics