

INVACARE®
JASMINE™
PREMIER SERIES FULL BODY LIFT



The **Invacare® Jasmine™ Full Body Lift** has optimum functionality and innovative features that provide safer patient handling within a 500 lb. weight capacity. Caregivers can rely on safer operation and patients on maximum comfort during everyday transfers. Beyond that, the Jasmine lift expertly functions in more challenging situations with a rotating 360° 6-point hanger bar and 16.4" minimum height range, essential for assisting patients who may have fallen to the floor.

The stable spring-loaded hanger bar produces an audible click when the sling is engaged. This assures the caregiver that the sling is properly attached. That means less stress for the caregiver, a stable transfer for the patient and less likelihood of injury for both. Optimized maneuverability is key for a comfortable transfer. Jasmine lifts feature a battery-powered boom and wide-opening base, so caregivers can easily walk behind them and wheelchairs can be accommodated up front. The push handle accommodates any caregiver, and with a lighter-weight aluminum and steel construction, control in movement is straightforward and trouble-free.

Features

Lifting Range: This lift expertly functions in more challenging situations with a rotating 360° 6-point hanger bar and 16.4" minimum height range, essential for assisting patients who may have fallen to the floor

Maneuverability: The push handle accommodates any caregiver. And with a lighter-weight aluminum and steel construction, control in movement is straightforward and trouble-free

Ease For Caregivers: Jasmine lifts feature a battery-powered boom and wide-opening base, both controlled through the hand pendant

Security: The stable spring-loaded hanger bar produces an audible tone when the sling is engaged. This assures the caregiver that the sling is properly attached

Cradle

XL Cradle

Scale : JLS5

Anti-Entrapment motor safety device



- 1. 360° 6 Point Hanger Bar
- 2. Hand Pendant
- 3. Push Handles
- 4. Cradle
- 5. Scale

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 11-243
111279

Model #	Height At Sling Hook Up	Base Dimensions	Casters	Weight	Weight Capacity	Low Battery Alarm	Emergency Controls
JASMINE	16.4"-74.6" H	4.8"H Clearance 46.7"W Open 25.8"W Closed 49"L	4"	98 lb.	500 lb.	Audio	Stop Up and Down
	Battery	Charge Time	Charger Input	Charger Output	Limited Warranty		
	24V DC Rechargeable 2 w/Unit	6 hr.	100-240V AC	29.5 DC	3 Years Lift 2 Years Electronics		

