

Nursing and care table MONA

The optimal solution for demanding care ...

The manifold demands of sheltered accommodation or homes for disabled people require highest demands on function and quality of care devices. In several hundreds accommodations therefore one trusts in our MONA table ...



Nursing and care table MONA

The nursing and care table MONA is an ideal support for people who require a high level of care. A large, adapted to the respective needs and requirements lying surface with high-quality padding is a basis for safe and comfortable care of a disabled person. The sturdy scissor-height-adjustment offers a necessary hold. Even for very anxious patients our system makes sure that neither the lying surface nor the table itself swings or yields. On request the nursing and care table MONA can be equipped with foldable side rails on each side. The side rails are padded to one's choice and even under high pressure do not fold offering an additional safety. A transfer of disabled persons can be implemented either by means of a mobile patient hoist - the table frame can be passed underneath - or directly and without difficulty from a wheelchair.

Product advantages nursing and care table MONA

- Mobile and easy to move
- Stable and sturdy chassis
- Stepless electrical or hydraulic height adjustment
- Easy to operate
- Applicable with all devices
- Different padding width and length possible
- Comfortable soft padding made out of artificial leather
- Highest safety against unintended rolling down
- Transfer in wheelchair height possible
- Max. patient weight 150 kg
- Optional: roll protection side bars for each side of the nursing and care table, padding for side rails, head part with catch adjustment