

INVACARE®  
**RELIANT™ 350**  
STAND-UP LIFT

---



---

The **Invacare® Reliant™ Series Stand-Up Lifts** are an ideal product for everyday patient handling, ensuring caregiver security and patient peace of mind during the transfer. Intended for use with weight-bearing patients as well as those needing rehabilitation support, the Invacare Reliant Stand-Up lift will be a valuable addition to your safety management system.

---

## Features

Innovative adjustment features adapt to a wide variety of body sizes and shapes

Low friction casters significantly improve rollability and contribute to caregiver safety

Multi-function slings allow for quick toileting, stand assist or full-support seated transfer

Adjustable leg support features a soft pad with lateral contour that secures legs in three positions for comfort at any height

Wide, non-slip footplate provides stable base for transfer

Under-bed clearance of 4.5" allows accessibility to most nursing home, hospital and homecare bed styles

Available with power opening base (Model # RPS350-2)

Removable battery eliminates downtime


Ergonomic features make it ideal for toilet transfers, as well as for dependent and rehabilitation purposes

Features include fast, simple sling attachment that protects against accidental disengagement; wide push handle; good access for clothing removal; and comfortable slings with extended back support that prevent pressure in the under-arm area

Adjustable leg strap option secures leg against knee pad for optimum positioning on footplate

Manufactured in an ISO13485 facility

## HCPCS Code: E0635



1. Electronics
2. Extra-large standing sling\* Model # R130
3. Extra-large transfer sling\* Model # R134

\*Complete instructions for lift usage and sling procedures are found in the Invacare operating manual. Only Invacare® slings and accessories should be used on Invacare® lifts.



**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012

Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.  
Form No. 11-248  
111284



Model #	Height At Sling Hook Up	Base Dimensions	Casters	Weight	Weight Capacity	Low Battery Alarm	Emergency Controls
RPS350-1	39.6"-63.7" H	4.5"H Clearance 36.8"W Open 25.7"W Closed 35.5"L	3" Front 5" Rear	115 lb.	350 lb.	Audio	Stop and Down
RPS350-2	39.6"-63.7" H	4.5"H Clearance 35.8"W Open 25.7"W Closed 35.5"L	3" Front 5" Rear	118 lb.	350 lb.	Audio	Stop and Down



	Battery	Charge Time	Charger Input	Charger Output	Limited Warranty
RPS350-1	N/A	N/A	N/A	N/A	3 Years Lift 1 Year Motor
RPS350-2	N/A	N/A	N/A	N/A	3 Years Lift 1 Year Motor

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.