

Diagnostic Solutions

# IQVITALS®



# Why Midmark Vitals?

Midmark IQvitals® is an affordable, Windows® based, touchscreen device that provides automated blood pressure, temperature, SpO<sub>2</sub> and pulse rate readings with a touch of the screen. IQvitals® can be used as a standalone device, connected to a computer using IQmanager® software or connect seamlessly to many leading EHR/EMRs. The intuitive user interface helps increase productivity, reduce errors and simplify staff training. Take your practice to the next level with Midmark IQvitals®.



## The path to better care begins with vitals.

Vitals are the beginning of the patient-caregiver interaction. They provide critical information related to changes in patient health and play an important role in a physician's treatment decisions. The integration of EMRs, new technologies and automated devices has a significant impact on the overall ease and efficiency of the process. So why are we still using the same vitals process we have been for the past 30 years?

**It's time for a change.**



Lightweight and portable with flexible mounting options to meet the needs of any clinical environment.



Automated capture of vital signs data improves workflow efficiency and helps reduce transcription errors.



Intuitive color touchscreen simplifies usability and reduces learning curves.

## IQvitals® technology by Midmark

Use as a standalone unit or connect to a computer using IQmanager® software or leading EHR/EMRs



All-in-one vital signs device (with options to add SpO<sub>2</sub> and digital scales) increases efficiency.



Seamless connectivity to the patient chart in the EMR or IQmanager® helps reduce errors.

# Intuitive, efficient user interface

The graphical user interface was designed by clinicians for clinicians. Information is displayed in an easy-to-read format with the most commonly used results most prominent for quick reference. The screen menus are simple to navigate and operate, allowing you to easily align with your settings to your practice's information flow.

BP results prominently displayed in the top corner

Temperature, SpO<sub>2</sub> and pulse rate framed with a soft background for quick identification

Prominent battery status to ensure uninterrupted care



Conveniently located BP presets

Clear indication of selected temperature mode

Automatically calculates BMI when height and weight is captured

## Control the vitals process from one easy screen

Operate the Midmark Barrier-Free® Exam Table with IQscale® and IQvitals® from one integrated control screen on the IQvitals® device, saving valuable time to make you more efficient.



Capture weight with one touch while the patient is comfortably seated on the exam table.




## Truly connected for the first time

Consider the future of healthcare—it will mean more connectivity, an increase in integrated data and the expansion of the EMR. And, in many ways, that future is already here. The good news is that innovations in equipment are now available to help you and your staff navigate this new technological terrain.

The Midmark 625 Barrier-Free® Examination Table with IQscale® and IQhub® allows you to connect IQvitals® and other devices directly to the table and port diagnostic data straight into the EMR. This makes gathering data and managing patient records easier, faster and more accurate. Plus, state-of-the-art connectivity gives you the ability to better meet mandatory Meaningful Use requirements.



## Features



# Featured Accessories

Giving you the tools you need to improve quality of care is so important to us.  
We want you to have everything you need to make the best use of your vitals device.



IQvitals® Mobile Cart



IQcart®



6231 Care Exchange® Workstation




IQvitals® Wall Mount



IQvitals® Equipment Pole Mount



IQvitals® Countertop Mount



Fairbanks® TeleWeigh™ Digital Scale



625 Table Mount for IQvitals® - IQhub®



625 Table Mount for IQvitals® - Non-IQhub®



Reusable Blood Pressure Cuffs



Turbo Temp® Probe Covers



Fairbanks® Scale Mount



Reusable SpO<sub>2</sub> Finger Sensor by Envitec®



IQscale® Serial Cable for IQvitals® - 6 foot



IQscale® Serial Cable for IQvitals® - 15 foot



IQvitals® Thermal Printer



IQvitals® Mobile Cart Printer Bracket

# Make your practice more efficient with an IQ Digital Diagnostic System

Workflow efficiency and improvements in patient care are just a few of the benefits that can result from implementing an IQ Digital Diagnostic System. Call us today and we'll put together a customized suite of solutions to help grow as your practice grows.

Diagnostic Products				
	EHR/EMR Connectivity	Single Platform (IQmanager®)	Standalone	Thin-client
IQecg®	■	■		■
IQholter®	■	■		■
IQholter® EX	■	■		■
IQholter® EP	■	■		■
IQspiro®	■	■		■
IQvitals®	■	■	■	■

## Minimum System Requirements

- Windows® 8, Windows® 7; Professional and Enterprise; 32-bit and 64-bit
- Intel® Core™ 2 Duo Processor E4300 or 64-bit processor or faster
- 2 GB of RAM
- 2 GB available hard disk space
- USB or Serial Port
- Mouse, keyboard and high resolution display monitor

*Windows is a registered trademark of Microsoft Corporation.*

*Intel Core is a trademark of Intel Corporation in the U.S. and/or other countries.*

Contact us or visit [midmark.com](http://midmark.com) for more detailed system requirements.

## Thin-client System Requirements

Please contact your local Midmark Sales Representative for detailed system requirements.

## Ordering Information

4-000-0500  
IQvitals® with Blood Pressure and Alaris® Turbo Temp® Thermometer

4-000-0510  
IQvitals® with Blood Pressure, Alaris® Turbo Temp® Thermometer and SpO<sub>2</sub>

*Alaris and Turbo Temp are registered trademarks of Cardinal Health, Inc.*

*Fairbanks and TeleWeigh are trademarks of Fairbanks Scales, Inc.*

## Accessories

IQvitals® Mobile Cart  
3-004-2000

IQcart®  
3-004-1000

6231 Care Exchange® Workstation  
6231-001-xxx

IQvitals® Wall Mount  
3-009-0003

IQvitals® Equipment Pole Mount  
3-004-2008

IQvitals® Countertop Mount  
3-009-0001

Fairbanks® TeleWeigh™ Digital Scale  
1-100-1603

Table Mount for IQvitals®  
9A551001 IQhub®

Table Mount for IQvitals®  
9A551002 Non-IQhub®

Reusable Blood Pressure Cuffs  
3-009-0068 Infant, 3-009-0070 Child, 3-009-0062 Small Adult, 3-009-0064 Adult, 3-009-0066 Large Adult, 3-009-0072 Adult Long, 3-009-0074 Large Adult Long 3-009-0076 Thigh

Turbo Temp® Probe Covers  
(20 covers per box) 3-009-0058

Fairbanks® Scale Mount  
3-004-2010

Reusable SpO<sub>2</sub> Finger Sensor by Envitec®  
3-009-0021 Small, 3-009-0020 Adult, 3-009-0026 SpO<sub>2</sub> 4' Extended Cable

IQscale® Serial Cable for IQvitals®  
9A478001 6 foot length, coiled

IQscale® Serial Cable for IQvitals®  
9A478002 15 foot length, straight

IQvitals® Thermal Printer  
1-100-1605

IQvitals® Mobile Cart Printer Bracket  
3-004-2004

**Midmark is an ISO 13485 and ISO 9001 Certified Company**  
For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-800-365-8631  
Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8392  
or visit our website at [midmark.com](http://midmark.com)

© 2004 Midmark Corporation  
Midmark Corporation, Dayton, OH.  
Products subject to improvement changes without notice  
Litho in U.S.A. 21-58-0001 Rev. C (7/16)

