



SPIRO-S
SpiroPerfect® Module
(Laptop not included)

PC-Based SpiroPerfect® Spirometer

See and compare accurate results from pulmonary function tests in seconds. Clearly evaluate pulmonary obstruction or restriction, provide motivating feedback for smoking cessation patients. Show progress, and help pediatric patients comply with test procedures.

- Fast, reliable pre-/postbronchodilator FVC, SVC, and MVV testing with automatic interpretation and reversibility analysis
- Instantly verify quality and variability of test performance and compare previous tests using the trending tool
- Motivate children with incentive graphic
- Reduce risk of cross contamination with single-use disposable transducers
- Easily monitor, review, store and print on your computer with Welch Allyn CardioPerfect® Workstation Software
- Interface with EHR or patient management systems for onsite or offsite review of key data

SPIRO-S SpiroPerfect Module, CPWS, with 3-liter calibration syringe (Includes: sensor, pressure tubing, flow transducers, electronic operator's manual, latest CPWS Software CD, 2 training posters and quick reference card)

SPIRO-NS SpiroPerfect Module, CPWS, without 3-liter calibration syringe (Includes: sensor, pressure tubing, flow transducers, electronic operator's manual, latest CPWS Software CD, 2 training posters and quick reference card)

SPIRO-S-RS232 SpiroPerfect Module, CPWS, with 3-liter calibration syringe. Includes: CPWS spiro Software module, sensor, pressure tubing, 8 flow transducers, electronic manual, 2 training posters, and quick reference card)


SPIRO-NS-RS232 SpiroPerfect Module, CPWS, w/o 3-liter calibration syringe. Includes: CPWS spiro Software module, sensor, pressure tubing, 8 flow transducers, electronic manual, 2 training posters, and quick reference card



703418
Disposable Flow
Transducers



703480
Calibration Syringe



100680
Spirometry Nose Clips



703415
Pressure Tubing

PC-Based SpiroPerfect Spirometer: Accessories

703419 Disposable Flow Transducers, 100 pk, CPWS, CP 200

703418 Disposable Flow Transducers, 25 pk, CPWS, CP 200

703337 Effort Acceptability Poster

71038-3000 Spirometry Reference Chart

703480 Calibration Syringe

100680 Spirometry Nose Clips, (25/pack)

703415 Pressure Tubing, CPWS, CP 200, (2 m)

Viewing Licenses

101524 Viewing License

101526 Viewing License